



# Cole-Parmer®

## CYBERSCAN PH / ION 510 BENCH METER

~~List Price : Rs. 49,999/-~~

**Special Promotional Price : Rs. 38,500/-**

+

**Epoxy Body pH Electrode**

**FREE**

### HIGH PERFORMANCE AT AN ECONOMICAL COST



- Multi-parameter measurement
- Ion measurement capabilities with direct readout in ppm by use of Ion Selective Electrode (with BNC connector)
- 5 point calibration with auto-buffer recognition
- 50 memory data sets
- Self diagnosis with user friendly message codes
- Automatic temperature compensation
- USA or NIST buffer standards

**What's included:**

CyberScan pH/Ion 510 meter + Temperature probe + Integral electrode holder + Power adapter (230V)

#### Specifications

ISO9001:2000  
CERTIFIED SUPPLIER



3 year warranty

Slope Control: 80 to 120%  
Operating temperature: 0 to 50°C

Range	pH	0.00 to 14.00 pH
	mV	±1999 mV
	Temperature	0 to 100 °C
	Concentration	0.01 to 1999 ppm
Resolution	pH	0.01 pH
	mV	0.1 mV within ±199.9 mV, 1 mV beyond ±199 mV
	Temperature	0.1 °C
	Concentration	0.01 / 0.1 / 1 ppm
Accuracy	pH	±0.01 pH
	mV	0.2 mV within ±199.9 mV, 2 mV beyond ±199 mV
	Temperature	±0.5 °C
	Concentration	±5% of reading
Output	±1999 mV, 1mV/mV unit via mini pin recorder connector	
Display	Dual LCD with mode annunciators	
Power	110/220 VAC	

**Terms & Conditions:**

Prices quoted are Ex-works warehouse  
Taxes & other levies extra as applicable  
Offer valid only upto 31st March 2010



**GO to ColeParmer.in**  
Complete the system by  
ordering ISE of your choice



**For more information contact**

**Cole-Parmer India**

(Div. of Scientific Products India Pvt Ltd)  
Godrej Coliseum 101 A - 101 B,  
Somaiya Hospital Road, Sion (E),  
Mumbai - 400022, India  
Tel: +91-22-67162209  
Fax: +91-22-67162211  
[shraddha.damle@coleparmer.in](mailto:shraddha.damle@coleparmer.in)